

The Board of Pharmacy has received notice of the following product recall:

Product: BD Vacutainer® Push Button Blood Collection Set with Pre-Attached Holder

Catalog Number: 368656

Lot Number	UDI (GTIN, DI + PI)	Exp. Date
9079767	(01)50382903686563 (17)210331 (10)9079767(30)0100	3/31/2021
9079770	(01)50382903686563 (17)210331 (10)9079770(30)0100	3/31/2021
9081829	(01)50382903686563 (17)210331 (10)9081829(30)0100	3/31/2021
9086811	(01)50382903686563 (17)210331 (10)9086811(30)0100	3/31/2021
9086812	(01)50382903686563 (17)210331 (10)9086812(30)0100	3/31/2021
9091607	(01)50382903686563 (17)210331 (10)9091607(30)0100	3/31/2021
9091626	(01)50382903686563 (17)210331 (10)9091626(30)0100	3/31/2021
9094659	(01)50382903686563 (17)210331 (10)9094659(30)0100	3/31/2021
9094661	(01)50382903686563 (17)210331 (10)9094661(30)0100	3/31/2021
9098547	(01)50382903686563 (17)210430 (10)9098547(30)0100	4/30/2021
9100711	(01)50382903686563 (17)210430 (10)9100711(30)0100	4/30/2021
9100712	(01)50382903686563 (17)210430 (10)9100712(30)0100	4/30/2021

BD is recalling for certain lots of BD Vacutainer® Push Button Blood Collection Set with Pre-Attached Holder, Catalog# 368656, 23G x 3/4" x 12" as listed above. BD has confirmed that these lots may exhibit separation of front and rear barrels upon activation of the safety feature that retracts the needle. Barrel separation may result in leakage of blood or an exposed needle, potentially resulting in exposure to blood borne pathogens.

The Board of Pharmacy strongly encourages pharmacies to immediately review their quality assurance and recall policies and procedures to determine if any corrective action is required.